



# RNAseq in Nextflow

Yuan Zhang  
HPC Facility  
October 12, 2023



# Topics

- Nextflow basics
- A simple nextflow example
- RNAseq in nextflow


# What is Nextflow?

It is a workflow management software which enables the writing of scalable and reproducible scientific workflows.

Scalability - It implements a Domain Specific Language (DSL) that simplifies the implementation and running of workflows on cloud or HPC infrastructure.


Reproducibility - It integrates various software packages and environment management system such as Docker, Singularity, and Conda.

# Nextflow Diagram



# Nextflow Basics - Process

- Describes a task to be run.
- The script can be written in any language (Bash, Perl, Python, R, etc.)
- One task for each input set.



## Example - Process

```
process NUM_LINES {  
    input: path read  
    output: stdout  
    script:  
    """  
    printf '${read}:'  
    gunzip -c ${read} | wc -l  
    """  
}
```

# Nextflow Basics - Channel

Stores the data as input or output of the processes. Processes communicate with each other through channels.

Value channel – can be used multiple times

- `Channel.value`

Queue channel – FIFO queues, used only once.

- `Channel.of`
- `Channel.fromList`
- `Channel.fromPath`
- `Channel.fromFilePairs`
- `Channel.fromSRA`

# Example - Channel

Channel contains one data set

```
groovy> input_ch =  
    Channel.fromPath("~/data/fastqs/SRR12005075.fastq.gz")  
groovy> input_ch.view()  
  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005075.fastq.gz
```


Channel contains multiple data set

```
groovy> input_ch =  
    Channel.fromPath("~/data/fastqs/*.fastq.gz")  
groovy> input_ch.view()  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005057.fastq.gz  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005058.fastq.gz  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005059.fastq.gz  
...  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005073.fastq.gz  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005074.fastq.gz  
/home/gdhpcgroup/yxz006/data/fastqs/SRR12005075.fastq.gz
```



# Nextflow Basics - Workflow

- Defines a sequence of tasks that processes a set of data
- The execution environment (local, slurm, aws, etc.) is defined in config file.
- Can be run on local computer, HPC cluster, or cloud.



## Example - Workflow

```
params.input =  
"~/data/fastqs/SRR12005075.fastq.gz"  
  
workflow {  
    input_ch = Channel.fromPath(params.input)  
    NUM_LINES(input_ch)  
    NUM_LINES.out.view()  
}
```

# A Simple Workflow – wc.nf

```
1 nextflow.enable.dsl=2
2
3 params.input = "~/data/fastqs/SRR12005075.fastq.gz"
4
5 workflow {
6     input_ch = Channel.fromPath(params.input)
7     NUM_LINES(input_ch)
8     NUM_LINES.out.view()
9 }
10
11 process NUM_LINES {
12     input: path read
13     output: stdout
14     script:
15         """
16         printf '${read}: '
17         gunzip -c ${read} | wc -l
18         """
19 }
```

# Run the Workflow on Local Computer

```
[yxz006@r1pl-hpcf-log02 wc]$ ml nextflow
```

```
[yxz006@r1pl-hpcf-log02 wc]$ nextflow run wc.nf
```

```
N E X T F L O W ~ version 22.10.6
```

```
Launching `wc.nf` [zen_lorenz] DSL2 - revision: 28933bf913
```

```
executor > local (1)
```

```
[b1/b70dc6] process > NUM_LINES... [100%] 1 of 1 ✓
```

```
SRR12005075.fastq.gz: 71499804
```

# Run the Workflow on HPC Cluster

```
[yxz006@r1pl-hpcf-log02 wc]$ cat nextflow.config  
process.executor = 'slurm'
```

```
[yxz006@r1pl-hpcf-log02 wc]$ nextflow run wc.nf  
N E X T F L O W ~ version 22.10.6  
Launching `wc.nf` [infallible_franklin] DSL2 - revision: 28933bf913  
executor > slurm (1)  
[35/46dd0f] process > NUM_LINES (1) [100%] 1 of 1 ✓  
SRR12005075.fastq.gz: 71499804
```

# Change the Input Parameter

```
[yxz006@r1pl-hpcf-log02 wc]$ nextflow run wc.nf -ansi-log false \  
    --input="/home/gdhpcgroup/yxz006/data/fastqs/*.gz"
```

```
...
```

```
Launching `wc.nf` [suspicious_roentgen] DSL2 - revision: 0e75a1e416
```

```
[a7/24860b] Submitted process > NUM_LINES (1)
```

```
...
```

```
[a8/25498a] Submitted process > NUM_LINES (23)
```

```
[b5/56c3de] Submitted process > NUM_LINES (24)
```

```
SRR12005057.fastq.gz: 71792404
```

```
...
```

```
SRR12005065.fastq.gz: 71545878
```

```
...
```

# The Resume Option

```
[yxz006@r1pl-hpcf-log02 wc]$ nextflow run wc.nf -ansi-log false \  
  --input="/home/gdhpcgroup/yxz006/data/fastqs/*.gz" \  
  -resume
```

...

```
Launching `wc.nf` [extravagant_marconi] DSL2 - revision: 64b7f7c6e0  
[a7/24860b] Cached process > NUM_LINES (1)
```


...

```
[a8/25498a] Cached process > NUM_LINES (23)  
[b5/56c3de] Submitted process > NUM_LINES (24)  
SRR12005057.fastq.gz: 71792404
```

...

```
SRR12005065.fastq.gz: 71545878  
SRR12005066.fastq.gz: 71833268
```

# Nextflow Process Flow Diagram





# A Simple Workflow – wc.nf

```
1 nextflow.enable.dsl=2
2
3 params.input = "~/data/fastqs/SRR12005075.fastq.gz"
4
5 workflow {
6     input_ch = Channel.fromPath(params.input)
7     NUM_LINES(input_ch)
8     NUM_LINES.out.view()
9 }
10
11 process NUM_LINES {
12     input: path read
13     output: stdout
14     script:
15         """
16         printf '${read}: '
17         gunzip -c ${read} | wc -l
18         """
19 }
```

# Matt's RNAseq Pipeline

- Input file
  - `misc/MiiceIiiiiinSpaaaaaaaaaaaaace.csv`
- Slurm scripts
  - `sbatchCmds/getSRADData.sh`
  - `sbatchCmds/fastqc.sh`
  - `sbatchCmds/align_hisat.sh`
  - `sbatchCmds/featureCounts.sh`
- Workflow
  - `MainRNAseqAnalysis.Rmd`

# Converting the Pipeline to Nextflow

- `misc/MiiceIiiinSpaaaaaaaaaaaaace.csv`
- `sbatchCmds/getSRADData.sh`
  - `getSRADData.nf`
- `sbatchCmds/fastqc.sh`
  - `fastqc.nf`
- `sbatchCmds/align_hisat.sh`
  - `align_hisat.nf`
- `sbatchCmds/featureCounts.sh`
  - `featureCounts.nf`
- `MainRNAseqAnalysis.Rmd`
  - `rnaseq.nf` + `MainRNAseqAnalysis.R` + `pca_de.nf`

# Nextflow Script - getSRAData.nf

```
process getSRAData {
    input:
    val (sra_number)

    output:
    path("*.gz"),          emit: fastq

    shell:
    '''
    module purge
    module load SRAToolkit/3.0.1
    prefetch -O ./ !{sra_number}

    fasterq-dump -e 4 -S -O ./ ./!{sra_number}/!{sra_number}.sra
    pigz ./!{sra_number}.fastq
    '''
}
```

# Nextflow Config – nextflow.config

```
process {  
    executor = 'slurm'  
    cpus = 10  
    time = '12hours'  
    queue = 'general'  
    memory = ''  
    email = ''  
  
    withName: getSRADData {  
        cpus = 4  
        time = '24hours'  
    }  
  
    withName: fastqc {  
        cpus = 4  
    }  
}
```

# Nextflow Script – rnaseq.nf

```
params.samples = "$params.samplerdir/MiiiceIiiiiinSpaaaaaaaaaaaaace.csv"
workflow RNASEQ {
    ch_samples = Channel.fromPath(params.samples)
    ch_samples.splitCsv ( header:true, sep:', ' )
    .map { it.Run }
    .set {sra_array}

    getSRADData (sra_array)

    ch_fastq = getSRADData.out.fastq
    fastqc ( ch_fastq )

    align_hisat (ch_fastq)

    ch_bam = align_hisat.out.bam.collect()
    featureCounts (ch_bam)

    ...
}
```

# Nextflow Script – main.nf

```
nextflow.enable.dsl = 2

// INPUT DIRECTORIES
params.projDir = "/home/gdhpcgroup/yxz006/training/2023/rnaseq/"
params.samples = "$params.projDir/misc/MiiiceIiiiiinSpaaaaaaaaaaaaace.csv"
params.fastqs = "$params.projDir/input/fastqs"

// OUTPUT DIRECTORIES
params.outdir = "$params.projDir/output"
params.fastqc = "$params.outdir/fastqc"
params.multiqc = "$params.outdir/multiqc"
params.align = "$params.outdir/align"
params.geneCounts = "$params.outdir/geneCounts"
params.figures = "$params.outdir/figures"
params.rdata = "$params.outdir/rdata"
params.de = "$params.outdir/de"
params.gsea = "$params.outdir/gsea"

// RUN WORKFLOW RNASEQ
include { RNASEQ } from './workflows/rnaseq'
workflow {
    RNASEQ ()
}
```

# Nextflow Script – pca\_de.nf

```
process pca_de {
    input:
        path ("*multiqc.log")

    shell:
        '''
        mkdir -p !{params.figures} !{params.rdata} !{params.de}
        !{params.gsea}

        ml purge
        ml GCC/9.3.0 OpenMPI/4.0.3 R/4.2.2

        workDir=$(pwd)
        cd !{params.projDir}
        Rscript ./MainRNAseqAnalysis.R

        # save the slurm output
        cd ${workDir}
        cat .command.log > ${SLURM_JOB_ID}_pca_de.log
        '''
}
```




# Run RNAseq in Nextflow

```
[yxz006@r1pl-hpcf-log02 rnaseq]$ nextflow run main.nf \  
    -resume \  
    -ansi-log false \  
    -with-dag flowchart.png \  
    -with-report resource_utilization.html
```

# RNAseq pipeline in DAG Format

DAG – Directed acyclic graph



# Resource Utilization Report

## Nextflow workflow report

[sleepy\_jepsen] (*resumed run*)

Workflow execution completed successfully!

### Run times

11-Oct-2023 11:27:58 - 11-Oct-2023 11:38:24 (duration: **10m 26s**)

75 succeeded

### Nextflow command

```
nextflow run -resume -ansi-log false main.nf -with-dag flowchart.png -with-report resource_utilization.html
```

|                   |  |
|-------------------|--|
| CPU-Hours         | 23.7   |
| Launch directory  | /gpfs0/home1/gdhpcgroup/yxz006/training/2023/rnaseq      |
| Work directory    | /gpfs0/home1/gdhpcgroup/yxz006/training/2023/rnaseq/work |
| Project directory | /gpfs0/home1/gdhpcgroup/yxz006/training/2023/rnaseq      |
| Script name       | main.nf  |
| Script ID         | 5751fd75e3e69822dc12af8f0d4f5389                         |
| Workflow session  | 0ca5aa22-7416-4ba5-ab73-0b2c7679eb05                     |
| Workflow profile  | standard   |
| Nextflow version  | version 22.10.6, build 5843 (23-01-2023 23:20 UTC)       |

# Resource Utilization Report - 2


## Resource Usage

These plots give an overview of the distribution of resource usage for each process.

### CPU

Raw Usage

% Allocated



# Resource Utilization Report - 3


## Memory

Physical (RAM)

Virtual (RAM + Disk swap)

% RAM Allocated

Physical Memory Usage



# Resource Utilization Report - 4

## Job Duration

Raw Usage

% Allocated



# Benefits of Running RNAseq in Nextflow

## Integrate various software packages

- QC, trimming, alignment, featurecounts, multiqc, differential expression analysis

## Reproducible using environment management system

- Docker, Singularity, Conda

## Portable

- local, Slurm, AWS

## Continuous checkpoints

- Run the pipeline from the last successfully executed steps.

# References:

- <https://training.seqera.io/>
- <https://carpentries-incubator.github.io/workflows-nextflow/>
- [https://github.com/MVesuviusC/bioinformatics\\_meeting](https://github.com/MVesuviusC/bioinformatics_meeting)